

产品名称: RCC1 Rabbit Monoclonal Antibody
产品货号: AMRe03184



产品概述 (Summary)

产品名称 (Production Name)	RCC1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RCC1
别名 (Alternative Names)	RCC1; RCC1-I; SNHG3-RCC1; CHC1
基因 ID (Gene ID)	1104
蛋白 ID (SwissProt ID)	P18754.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

Guanine-nucleotide releasing factor that promotes the exchange of Ran-bound GDP by GTP. Involved in the regulation of


产品名称: RCC1 Rabbit Monoclonal Antibody
产品货号: AMRe03184




onset of chromosome condensation in the S phase. Binds both to the nucleosomes and double-stranded DNA.

研究领域 (Research Area)

图片 (Image Data)



Immunocytochemistry analysis of RCC1 (green) in HeLa using RCC1 antibody, and DAPI (blue).



Western blot analysis of RCC1 in HeLa, MCF-7, HepG2, U251, HL-60 lysates using RCC1 antibody.

注意事项 (Note)

For research use only .