

产品名称: RhoGDI Rabbit Monoclonal Antibody
产品货号: AMRe03186



产品概述 (Summary)

产品名称 (Production Name)	RhoGDI Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARHGDI A
别名 (Alternative Names)	ARHGDI A; GDIA1; Rho GDP-dissociation inhibitor 1; Rho GDI 1; Rho-GDI alpha
基因 ID (Gene ID)	396
蛋白 ID (SwissProt ID)	P52565.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 26 kDa

研究背景 (Background)

Regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them, and the

产品名称: RhoGDI Rabbit Monoclonal Antibody
产品货号: AMRe03186




subsequent binding of GTP to them.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of RhoGDI (green) in Hela using RhoGDI antibody, and DAPI (blue).



Western blot analysis of RhoGDI in 3T3, mouse brain lysates using RhoGDI antibody.



Western blot analysis of RhoGDI in Hela, A549, HL-60, C6 lysates using RhoGDI antibody

产品名称: RhoGDI Rabbit Monoclonal Antibody
产品货号: AMRe03186



注意事项 (Note)

For research use only .