

产品名称: RNF14 Rabbit Monoclonal Antibody
产品货号: AMRe03187



产品概述 (Summary)

产品名称 (Production Name)	RNF14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF14
别名 (Alternative Names)	RNF14; ARA54; HRIHFB2038; E3 ubiquitin-protein ligase RNF14; Androgen receptor-associated protein 54; HFB30; RING finger protein 14; Triad2 protein
基因 ID (Gene ID)	9604
蛋白 ID (SwissProt ID)	Q9UBS8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

研究背景 (Background)


产品名称: RNF14 Rabbit Monoclonal Antibody
产品货号: AMRe03187




Might act as an E3 ubiquitin-protein ligase which accepts ubiquitin from specific E2 ubiquitin-conjugating enzymes and then transfers it to substrates, which could be nuclear proteins.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RNF14 in mouse brain, rat brain lysates using RNF14 antibody.



Western blot analysis of RNF14 in A549, HL-60, U2OS lysates using RNF14 antibody.

注意事项 (Note)

For research use only .