

产品名称: SFT Rabbit Monoclonal Antibody
产品货号: AMRe03202



产品概述 (Summary)

产品名称 (Production Name)	SFT Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE2D1
别名 (Alternative Names)	ubiquitin conjugating enzyme E2 D1; SFT; UBC5; UBC4/5; UBC5A; E2(17)KB1
基因 ID (Gene ID)	7321
蛋白 ID (SwissProt ID)	P51668.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)


产品名称: SFT Rabbit Monoclonal Antibody
产品货号: AMRe03202




Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of SFT in Raw264.7 lysates using SFT antibody.



Western blot analysis of SFT in HeLa, HL-60 lysates using SFT antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SFT antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: SFT Rabbit Monoclonal Antibody
产品货号: AMRe03202



注意事项 (Note)

For research use only .