

产品名称: SIRP alpha Rabbit Monoclonal Antibody
产品货号: AMRe03203



产品概述 (Summary)

产品名称 (Production Name)	SIRP alpha Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SIRPA SIRPA; BIT; MFR; MYD1; PTPNS1; SHPS1; SIRP; Tyrosine-protein phosphatase non-receptor type substrate 1; SHP substrate 1; SHPS-1; Brain Ig-like molecule with tyrosine-based activation motifs; Bit; CD172 antigen-like family member A; Inhibito
别名 (Alternative Names)	
基因 ID (Gene ID)	140885
蛋白 ID (SwissProt ID)	P78324.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 85 kDa

产品名称: SIRP alpha Rabbit Monoclonal Antibody
产品货号: AMRe03203




研究背景 (Background)

Immunoglobulin-like cell surface receptor for CD47. Acts as docking protein and induces translocation of PTPN6, PTPN11 and other binding partners from the cytosol to the plasma membrane. Supports adhesion of cerebellar neurons, neurite outgrowth and glial cell attachment.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of SIRP alpha in A549, MCF-7, 293T lysates using SIRP alpha antibody.



Immunocytochemistry analysis of SIRP alpha (green) in HeLa using SIRP alpha antibody, and DAPI (blue)

注意事项 (Note)

For research use only .