

产品名称: SKP2 Rabbit Monoclonal Antibody
产品货号: AMRe03204



产品概述 (Summary)

产品名称 (Production Name)	SKP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SKP2
别名 (Alternative Names)	SKP2; FBXL1; S-phase kinase-associated protein 2; Cyclin-A/CDK2-associated protein p45; F-box protein Skp2; F-box/LRR-repeat protein 1; p45skp2
基因 ID (Gene ID)	6502
蛋白 ID (SwissProt ID)	Q13309.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 48 kDa

研究背景 (Background)

产品名称: SKP2 Rabbit Monoclonal Antibody
产品货号: AMRe03204




Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunocytochemistry analysis of SKP2 (green) in hela using SKP2 Rabbit mAb, and DAPI(blue).

注意事项 (Note)

For research use only .