

产品名称: SMARCD1 Rabbit Monoclonal Antibody
产品货号: AMRe03207



产品概述 (Summary)

产品名称 (Production Name)	SMARCD1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMARCD1
别名 (Alternative Names)	CSS11; Rsc6p; BAF60A; CRACD1
基因 ID (Gene ID)	6602
蛋白 ID (SwissProt ID)	Q96GM5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 58 kDa

研究背景 (Background)

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-

产品名称: SMARCD1 Rabbit Monoclonal Antibody
产品货号: AMRe03207




nucleosome topology).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of SMARCD1 in HeLa, A549, HL-60 lysates using SMARCD1 antibody.

注意事项 (Note)

For research use only .