

产品名称: SPT5 Rabbit Monoclonal Antibody
产品货号: AMRe03211



产品概述 (Summary)

产品名称 (Production Name)	SPT5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SUPT5H
别名 (Alternative Names)	hSPT5; SPT5; SPT5H; supt5h; Tat CT1
基因 ID (Gene ID)	6829
蛋白 ID (SwissProt ID)	O00267.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 121 kDa; Observed MW: 150 kDa

研究背景 (Background)

Component of the DRB sensitivity-inducing factor complex (DSIF complex), which regulates mRNA processing and

产品名称: SPT5 Rabbit Monoclonal Antibody
产品货号: AMRe03211




transcription elongation by RNA polymerase II.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Immunocytochemistry analysis of SPT5 (green) in HeLa using SPT5 antibody, and DAPI (blue).



Western blot analysis of SPT5 in HeLa, A549, HL-60 lysates using SPT5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using SPT5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: SPT5 Rabbit Monoclonal Antibody
产品货号: AMRe03211



注意事项 (Note)

For research use only .