

产品名称: STAT6 Rabbit Monoclonal Antibody
产品货号: AMRe03215



产品概述 (Summary)

产品名称 (Production Name)	STAT6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAT6
别名 (Alternative Names)	STAT6; Signal transducer and activator of transcription 6; IL-4 Stat
基因 ID (Gene ID)	6778
蛋白 ID (SwissProt ID)	P42226.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 110 kDa

研究背景 (Background)

STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the

产品名称: STAT6 Rabbit Monoclonal Antibody
产品货号: AMRe03215




expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of STAT6 in HeLa, U2OS lysates using STAT6 antibody.

注意事项 (Note)

For research use only .