

产品名称: TACC3 Rabbit Monoclonal Antibody
产品货号: AMRe03221



产品概述 (Summary)

产品名称 (Production Name)	TACC3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TACC3
别名 (Alternative Names)	TACC3; ERIC1; Transforming acidic coiled-coil-containing protein 3; ERIC-1
基因 ID (Gene ID)	10460
蛋白 ID (SwissProt ID)	Q9Y6A5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 90 kDa; Observed MW: 140 kDa

研究背景 (Background)

TACC family members, TACC1, TACC2, and TACC3, map very closely to the corresponding FGFR1, FGFR2, FGFR3 genes on

产品名称: TACC3 Rabbit Monoclonal Antibody
产品货号: AMRe03221




chromosomes 4,8, and 10. Subsequently, since they are phylogenetically related, it is proposed that TACC and FGFR have similar roles in cell growth and differentiation.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of TACC3 in A549, HL-60, U251, U87-MG lysates using TACC3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TACC3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .