

**产品名称: Tissue Factor Pathway Inhibitor Rabbit
Monoclonal Antibody**
产品货号: AMRe03229



产品概述 (Summary)

产品名称 (Production Name)	Tissue Factor Pathway Inhibitor Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TFPI
别名 (Alternative Names)	EPI; TFI; LACI; TFPI1
基因 ID (Gene ID)	7035
蛋白 ID (SwissProt ID)	P10646.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 35 kDa; Observed MW: 40-55 kDa

研究背景 (Background)

**产品名称: Tissue Factor Pathway Inhibitor Rabbit
Monoclonal Antibody**
产品货号: AMRe03229




Inhibits factor X (X(a)) directly and, in a Xa-dependent way, inhibits VIIa/tissue factor activity, presumably by forming a quaternary Xa/LACI/VIIa/TF complex. It possesses an antithrombotic action and also the ability to associate with lipoproteins in plasma.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of TFPI in A549, HL-60, U2OS lysates using TFPI antibody.

注意事项 (Note)

For research use only .