

产品名称: Thrombospondin 1 Rabbit Monoclonal Antibody
产品货号: AMRe03231



产品概述 (Summary)

产品名称 (Production Name)	Thrombospondin 1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THBS1
别名 (Alternative Names)	TSP; THBS; TSP1; TSP-1; THBS-1
基因 ID (Gene ID)	7057
蛋白 ID (SwissProt ID)	P07996.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 129 kDa; Observed MW: 170 kDa

研究背景 (Background)

The protein encoded by this gene is a subunit of a disulfide-linked homotrimeric protein. This protein is an adhesive

产品名称: Thrombospondin 1 Rabbit Monoclonal Antibody
产品货号: AMRe03231




glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein can bind to fibrinogen, fibronectin, laminin, type V collagen and integrins alpha-V/beta-1. This protein has been shown to play roles in platelet aggregation, angiogenesis, and tumorigenesis. [provided by RefSeq, Jul 2008]


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Thrombospondin 1 in MCF-7, 293T, C2C12 lysates using Thrombospondin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Thrombospondin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .