

**产品名称: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody**
产品货号: AMRe03250



产品概述 (Summary)

产品名称 (Production Name)	Ubiquitin Conjugation Factor E4 B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBE4B
别名 (Alternative Names)	E4; UFD2; HDNB1; UBOX3; UFD2A
基因 ID (Gene ID)	10277
蛋白 ID (SwissProt ID)	O95155.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 146 kDa; Observed MW: 146 kDa

研究背景 (Background)


**产品名称: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody**
产品货号: AMRe03250




Ubiquitin-protein ligase that probably functions as an E3 ligase in conjunction with specific E1 and E2 ligases (By similarity).

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of UBE4B in PC-12, C6 lysates using Ubiquitin Conjugation Factor E4 B antibody.



Western blot analysis of UBE4B in HeLa, A549, HL-60, U2OS, U897MG lysates using UBE4B antibody.



Immunocytochemistry analysis of UBE4B (green) in MCF-7 using UBE4B antibody, and DAPI (blue)

**产品名称: Ubiquitin Conjugation Factor E4 B Rabbit
Monoclonal Antibody**
产品货号: AMRe03250



注意事项 (Note)

For research use only .