

产品名称: USP22 Rabbit Monoclonal Antibody
产品货号: AMRe03253



产品概述 (Summary)

产品名称 (Production Name)	USP22 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP22
别名 (Alternative Names)	Ubiquitin carboxyl terminal hydrolase 22; Usp22; USP3L
基因 ID (Gene ID)	23326
蛋白 ID (SwissProt ID)	Q9UPT9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa; Observed MW: 60 kDa

研究背景 (Background)

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the

产品名称: USP22 Rabbit Monoclonal Antibody
产品货号: AMRe03253




deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of USP22 in F9 lysates using USP22 antibody.



Western blot analysis of USP22 in A549, HL-60, U2OS, C6 lysates using USP22 antibody.

注意事项 (Note)

For research use only .