

产品名称: DAP kinase 3 Rabbit Monoclonal Antibody
产品货号: AMRe03261



产品概述 (Summary)

产品名称 (Production Name)	DAP kinase 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DAPK3
别名 (Alternative Names)	DAPK3; ZIPK; Death-associated protein kinase 3; DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase; Zipper-interacting protein kinase; ZIP-kinase
基因 ID (Gene ID)	1613
蛋白 ID (SwissProt ID)	O43293.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa; Observed MW: 53 kDa

研究背景 (Background)

产品名称: DAP kinase 3 Rabbit Monoclonal Antibody
产品货号: AMRe03261




Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ZIP Kinase in HeLa, A549, MCF-7, U251 lysates using DAP kinase 3 antibody.

注意事项 (Note)

For research use only .