

产品名称: Acetyl-Histone H4 (Lys5) Rabbit Monoclonal Antibody
产品货号: AMRe03268



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Histone H4 (Lys5) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H4C1
别名 (Alternative Names)	H4K5ac; H4F2; HIST; HIST1H4A; Histone H4.AC-H4K5; H4K5
基因 ID (Gene ID)	121504
蛋白 ID (SwissProt ID)	P62805.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 11 kDa; Observed MW: 11 kDa

研究背景 (Background)

产品名称: Acetyl-Histone H4 (Lys5) Rabbit Monoclonal Antibody
产品货号: AMRe03268




The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a member of the histone H4 family.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Immunocytochemistry analysis of Acetyl-Histone H4 (Lys5) (green) in HeLa using Acetyl-Histone H4 (Lys5) antibody, and DAPI (blue).




Western blot analysis of Acetyl-Histone H4 (Lys5) in 3T3 lysates using Acetyl-Histone H4 (Lys5) antibody. Observed band size: 11 kDa.

产品名称: Acetyl-Histone H4 (Lys5) Rabbit Monoclonal Antibody
产品货号: AMRe03268



Western blot analysis of Acetyl-Histone H4 (Lys5) in A549, HL-60, U2OS, U87-MG lysates using Acetyl-Histone H4 (Lys5) antibody



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Acetyl-Histone H4 (Lys5) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .