

产品名称: Acetyl-Histone H3 (Lys14) Rabbit Monoclonal Antibody
产品货号: AMRe03273



产品概述 (Summary)

产品名称 (Production Name)	Acetyl-Histone H3 (Lys14) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	H3C1
别名 (Alternative Names)	H3K14ac; H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J
基因 ID (Gene ID)	8350
蛋白 ID (SwissProt ID)	P68431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

产品名称: Acetyl-Histone H3 (Lys14) Rabbit Monoclonal Antibody
产品货号: AMRe03273




研究背景 (Background)

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Immunocytochemistry analysis of Acetyl-Histone H3 (Lys14) (green) in HeLa using Acetyl-Histone H3 (Lys14) antibody, and DAPI (blue).



Western blot analysis of Acetyl-Histone H3 (Lys14) in PC-12 lysates using Acetyl-Histone H3 (Lys14) antibody. Observed band size: 15 kDa.

产品名称: Acetyl-Histone H3 (Lys14) Rabbit Monoclonal Antibody
产品货号: AMRe03273



Western blot analysis of Acetyl-Histone H3 (Lys14) in A549, HL-60 lysates using Acetyl-Histone H3 (Lys14) antibody

注意事项 (Note)

For research use only .