

产品名称: Cyclin A1/A2 Rabbit Monoclonal Antibody
产品货号: AMRe03280



产品概述 (Summary)

产品名称 (Production Name)	Cyclin A1/A2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCNA2/CCNA1
别名 (Alternative Names)	CCN1; CCNA; CCNA1; CCNA2; CT146; Cyclin-A; Cyclin-A1; Cyclin-A2
基因 ID (Gene ID)	890/8900
蛋白 ID (SwissProt ID)	P20248/P78396.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 52,49 kDa; Observed MW: 52,49 kDa

研究背景 (Background)

May primarily function in the control of the germline meiotic cell cycle and additionally in the control of mitotic cell cycle in

产品名称: Cyclin A1/A2 Rabbit Monoclonal Antibody
产品货号: AMRe03280




some somatic cells. Essential for the control of the cell cycle at the G1/S (start) and the G2/M (mitosis) transitions.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Cyclin A1/A2 in A549, HL-60, U2OS lysates using Cyclin A1/A2 antibody.

注意事项 (Note)

For research use only .