

产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe03282



产品概述 (Summary)

产品名称 (Production Name)	PRPF31 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRPF31
别名 (Alternative Names)	Protein 61K; hPrp31; U4/U6 snRNP 61 kDa protein
基因 ID (Gene ID)	26121
蛋白 ID (SwissProt ID)	Q8WWY3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

研究背景 (Background)

Involved in pre-mRNA splicing as component of the spliceosome.


产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe03282




研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of PRPF31 in PC-12 lysates using PRPF31 antibody.



Western blot analysis of PRPF31 in C2C12, PC-3 lysates using PRPF31 antibody.



Immunocytochemistry analysis of PRPF31 (green) in HeLa using PRPF31 antibody, and DAPI (blue)

注意事项 (Note)

产品名称: PRPF31 Rabbit Monoclonal Antibody
产品货号: AMRe03282



For research use only .