

产品名称: Lysozyme Rabbit Monoclonal Antibody
产品货号: AMRe03289



产品概述 (Summary)

产品名称 (Production Name)	Lysozyme Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LYZ
别名 (Alternative Names)	LZM; LYZF1
基因 ID (Gene ID)	4069
蛋白 ID (SwissProt ID)	P61626.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)

Lysozymes have primarily a bacteriolytic function; those in tissues and body fluids are associated with the monocyte-

产品名称: Lysozyme Rabbit Monoclonal Antibody
产品货号: AMRe03289




macrophage system and enhance the activity of immunoagents.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Immunocytochemistry analysis of Lysozyme (green) in HEPG2 using Lysozyme antibody, and DAPI (blue).



Western blot analysis of Lysozyme in HL-60 lysates using Lysozyme antibody.



Western blot analysis of Lysozyme in Jurkat, C6, 3T3 lysates using Lysozyme antibody

产品名称: Lysozyme Rabbit Monoclonal Antibody
产品货号: AMRe03289



注意事项 (Note)

For research use only .