

产品名称: PLK1 Rabbit Monoclonal Antibody
产品货号: AMRe03292



产品概述 (Summary)

产品名称 (Production Name)	PLK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PLK1
别名 (Alternative Names)	PLK1; PLK; Serine/threonine-protein kinase PLK1; Polo-like kinase 1; PLK-1; Serine/threonine-protein kinase 13; STPK13
基因 ID (Gene ID)	5347
蛋白 ID (SwissProt ID)	P53350.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 68 kDa; Observed MW: 68 kDa

研究背景 (Background)

产品名称: PLK1 Rabbit Monoclonal Antibody
产品货号: AMRe03292




Required for recovery after DNA damage checkpoint and entry into mitosis. Required for kinetochore localization of BUB1B. Phosphorylates SGOL1. Required for spindle pole localization of isoform 3 of SGOL1 and plays a role in regulating its centriole cohesion function. Phosphorylates BORA, and thereby promotes the degradation of BORA. Contributes to the regulation of AURKA function. Regulates TP53 stability through phosphorylation of TOPORS.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of PLK1 in Jurkat, C6 lysates using PLK1 antibody.



Immunocytochemistry analysis of PLK1 (green) in HeLa using PLK1 antibody, and DAPI (blue)

注意事项 (Note)

For research use only .