

产品名称: CRCP Rabbit Monoclonal Antibody
产品货号: AMRe03297



产品概述 (Summary)

产品名称 (Production Name)	CRCP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRCP
别名 (Alternative Names)	C17; RCP; RCP9; RPC9; POLR3I; POLR3J; CGRPRCP; CGRP-RCP
基因 ID (Gene ID)	27297
蛋白 ID (SwissProt ID)	O75575.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 17 kDa; Observed MW: 17 kDa

研究背景 (Background)

DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates

产品名称: CRCP Rabbit Monoclonal Antibody
产品货号: AMRe03297




as substrates.


研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Immunocytochemistry analysis of CRCP (green) in HeLa using CRCP antibody, and DAPI (blue).




Western blot analysis of CRCP in 3T3 lysates using CRCP antibody.



Western blot analysis of CRCP in HeLa, C6 lysates using CRCP antibody

产品名称: CRCP Rabbit Monoclonal Antibody
产品货号: AMRe03297



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CRCP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .