

产品名称: Rabbit IgG (5B2) Mouse Monoclonal Antibody
产品货号: AMRe03532



产品概述 (Summary)

产品名称 (Production Name)	Rabbit IgG (5B2) Mouse Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Rabbit

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Rabbit IgG
别名 (Alternative Names)	Antibody; Ig gamma chain C region; Ig kappa-b4 chain C region; Ig kappa-b5 chain C region; Ig kappa-b9 chain C region; Ig lambda chain C region; K-BAS.
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	-.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 55,25 kDa; Observed MW: 55,25 kDa

研究背景 (Background)

产品名称: Rabbit IgG (5B2) Mouse Monoclonal Antibody
产品货号: AMRe03532




Constant region of immunoglobulin heavy chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of MEK2 in THP-1, Jurkat, Raji, 3T3 and Hela lysates using MEK2 antibody, then incubated with 201069(Rabbit IgG mAb).Anti-mouse IgG, HRPLinked antibody was used for detecting.



Western blot analysis of Rabbit IgG (5B2) in decreasing concentrations of total rabbit IgG reduced and denatured in 2X SDS loading buffer with DTT using IgG (HeavyChain Specific)antibody.



产品名称: Rabbit IgG (5B2) Mouse Monoclonal Antibody
产品货号: AMRe03532



Immunoprecipitation analysis of Rabbit IgG (5B2) in CHO-K1 lysates transfected with HAtagged NGAL using HATag antibody (lane 2). Lane 1 contains rabbit IgG IP control lysates. All lanes are probed with HATag antibody as the primary antibody. Panel A incubated with secondary antibodies anti-rabbit whole IgG, HRPlinked antibody for detecting. Panel B incubated with anti-rabbit IgG (HeavyChain Specific) mAb and the Anti-mouse whole IgG, HRPlinked antibody was used for detecting.

注意事项 (Note)

For research use only .