

产品名称: Cytochrome C Rabbit Monoclonal Antibody
产品货号: AMRe03749



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome C Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYCS
别名 (Alternative Names)	CYCS; CYC; Cytochrome c
基因 ID (Gene ID)	54205
蛋白 ID (SwissProt ID)	P99999.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 12 kDa; Observed MW: 12 kDa

研究背景 (Background)

CYCS Electron carrier protein. The oxidized form of the cytochrome c heme group can accept an electron from the heme

产品名称: Cytochrome C Rabbit Monoclonal Antibody
产品货号: AMRe03749




group of the cytochrome c1 subunit of cytochrome reductase. Cytochrome c then transfers this electron to the cytochrome oxidase complex, the final protein carrier in the mitochondrial electron-transport chain.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of Cytochrome C in mouse brain lysates using Cytochrome C antibody.



Western blot analysis of Cytochrome C in Jurkat, rat Brain, C6 lysates using Cytochrome C antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Cytochrome C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Cytochrome C Rabbit Monoclonal Antibody
产品货号: AMRe03749



注意事项 (Note)

For research use only .