

产品名称: Phospho-TrkA/B (Tyr490/Tyr516) Rabbit Monoclonal Antibody
产品货号: AMRe03752



产品概述 (Summary)

产品名称 (Production Name)	Phospho-TrkA/B (Tyr490/Tyr516) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NTRK2
别名 (Alternative Names)	NTRK2; TRKB; BDNF/NT-3 growth factors receptor; GP145-TrkB; Trk-B; Neurotrophic tyrosine kinase receptor type 2; TrkB tyrosine kinase; Tropomyosin-related kinase B
基因 ID (Gene ID)	4915
蛋白 ID (SwissProt ID)	Q16620.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa; Observed MW: 140 kDa

产品名称: Phospho-TrkA/B (Tyr490/Tyr516) Rabbit
Monoclonal Antibody
产品货号: AMRe03752




研究背景 (Background)

The family of Trk receptor tyrosine kinases consists of TrkA, TrkB and TrkC. While the sequence of these family members is highly conserved, they are activated by different neurotrophins: TrkA by NGF, TrkB by BDNF or NT4, and TrkC by NT3. TrkA regulates proliferation and is important for development and maturation of the nervous system. Point mutations, deletions and chromosomal rearrangements (chimeras) cause ligand-independent receptor dimerization and activation of TrkA.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Phospho-TrkA (tyr490)/TrkB (Tyr516) in K562, C6, 3T3 lysates using Phospho-TrkA/B (Tyr490/Tyr516) antibody.

注意事项 (Note)

For research use only .