

产品名称: CTGF Rabbit Monoclonal Antibody
产品货号: AMRe03767



产品概述 (Summary)

产品名称 (Production Name)	CTGF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCN2
别名 (Alternative Names)	CCN2; NOV2; HCS24; IGFBP8
基因 ID (Gene ID)	1490
蛋白 ID (SwissProt ID)	P29279.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 35 kDa

研究背景 (Background)

Major connective tissue chemoattractant secreted by vascular endothelial cells. Promotes proliferation and differentiation of

产品名称: CTGF Rabbit Monoclonal Antibody
产品货号: AMRe03767




chondrocytes. Mediates heparin- and divalent cation-dependent cell adhesion in many cell types including fibroblasts, myofibroblasts, endothelial and epithelial cells. Enhances fibroblast growth factor-induced DNA synthesis.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of CTGF in C6, HeLa lysates using CTGF antibody.

注意事项 (Note)

For research use only .