

产品名称: SQSTM1/p62 Rabbit Monoclonal Antibody
产品货号: AMRe03796



产品概述 (Summary)

产品名称 (Production Name)	SQSTM1/p62 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SQSTM1
别名 (Alternative Names)	p60; p62; A170; DMRV; OSIL; PDB3; ZIP3; p62B; NADGP; FTDALS3
基因 ID (Gene ID)	8878
蛋白 ID (SwissProt ID)	Q13501.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 62 kDa

研究背景 (Background)

Autophagy receptor that interacts directly with both the cargo to become degraded and an autophagy modifier of the

产品名称: SQSTM1/p62 Rabbit Monoclonal Antibody
产品货号: AMRe03796




MAP1 LC3 family. Required both for the formation and autophagic degradation of polyubiquitin-containing bodies, called ALIS (aggresome-like induced structures) and links ALIS to the autophagic machinery. Involved in midbody ring degradation. May regulate the activation of NFKB1 by TNF-alpha, nerve growth factor (NGF) and interleukin-1.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Immunocytochemistry analysis of SQSTM1/p62 (green) in U87-MG using SQSTM1/p62 antibody, and DAPI (blue).



Western blot analysis of SQSTM1 in K562, 3T3, HeLa lysates using SQSTM1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using SQSTM1 antibody. High-pressure

产品名称: SQSTM1/p62 Rabbit Monoclonal Antibody
产品货号: AMRe03796



and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .