

产品名称: Trk B Rabbit Monoclonal Antibody
产品货号: AMRe03802



产品概述 (Summary)

产品名称 (Production Name)	Trk B Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ntrk2
别名 (Alternative Names)	Tkrb; trkB; trk-B; GP145-TrkB/GP95-TrkB
基因 ID (Gene ID)	18212.0
蛋白 ID (SwissProt ID)	P15209.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa; Observed MW: 90-140 kDa

研究背景 (Background)

Receptor tyrosine kinase involved in the development and the maturation of the central and the peripheral nervous

产品名称: Trk B Rabbit Monoclonal Antibody
产品货号: AMRe03802




systems through regulation of neuron survival, proliferation, migration, differentiation, and synapse formation and plasticity. Receptor for BDNF/brain-derived neurotrophic factor and NTF4/neurotrophin-4. Alternatively can also bind NTF3/neurotrophin-3 which is less efficient in activating the receptor but regulates neuron survival through NTRK2.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of TrkB in mouse brain, rat brain lysates using Trk B antibody.

注意事项 (Note)

For research use only .