

产品名称: SLC31A1 Rabbit Monoclonal Antibody
产品货号: AMRe03812



产品概述 (Summary)

产品名称 (Production Name)	SLC31A1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat,Hamster

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC31A1
别名 (Alternative Names)	SLC31A1; COPT1; CTR1; solute carrier family 31 member 1
基因 ID (Gene ID)	1317
蛋白 ID (SwissProt ID)	O15431.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa; Observed MW: 26-34 kDa

研究背景 (Background)

The protein encoded by this gene is a high-affinity copper transporter found in the cell membrane. The encoded protein

产品名称: SLC31A1 Rabbit Monoclonal Antibody
产品货号: AMRe03812




functions as a homotrimer to effect the uptake of dietary copper.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of SLC31A1/CTR1 in mouse intestines lysates using SLC31A1 antibody.



Western blot analysis of SLC31A1/CTR1 in Hela, CHO-K1, rat Brain, Ramos lysates using SLC31A1/CTR1 antibody.

注意事项 (Note)

For research use only .