

产品名称: Phospho-alpha Tubulin (Tyr272) Rabbit Monoclonal Antibody
产品货号: AMRe03889



产品概述 (Summary)

产品名称 (Production Name)	Phospho-alpha Tubulin (Tyr272) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUBA1B
别名 (Alternative Names)	Tubulin alpha-1B chain; Alpha-tubulin ubiquitous; Tubulin K-alpha-1; Tubulin alpha-ubiquitous chain
基因 ID (Gene ID)	10376
蛋白 ID (SwissProt ID)	P68363.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa; Observed MW: 54 kDa

**产品名称: Phospho-alpha Tubulin (Tyr272) Rabbit
Monoclonal Antibody**
产品货号: AMRe03889




研究背景 (Background)

TUBA1A Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis on Jurkat lysates treated with pervanadate using Phospho-alpha Tubulin (Tyr272) antibody.

注意事项 (Note)

For research use only .