

产品名称: IRF5 Rabbit Monoclonal Antibody
产品货号: AMRe03898



产品概述 (Summary)

产品名称 (Production Name)	IRF5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IRF5
别名 (Alternative Names)	Interferon regulatory factor 5; Interferon regulatory factor 5 bone marrow variant; IRF 5; SLEB10
基因 ID (Gene ID)	3663
蛋白 ID (SwissProt ID)	Q13568.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 56 kDa

研究背景 (Background)

产品名称: IRF5 Rabbit Monoclonal Antibody
产品货号: AMRe03898




Play an important role in pathogen defense, autoimmunity, lymphocyte development, cell growth, and susceptibility to transformation. Studies have shown that genetic variants of IRF-5 have been associated with disorders where the IFN pathway is abnormally activated, such as systemic lupus erythematosus.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of IRF5 in THP-1 lysates using IRF5 antibody.

注意事项 (Note)

For research use only .