

产品名称: RIP Rabbit Monoclonal Antibody
产品货号: AMRe03918



产品概述 (Summary)

产品名称 (Production Name)	RIP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ripk1
别名 (Alternative Names)	Ripk1, Rinp, Rip, Rinp, Rip, Receptor-interacting serine/threonine-protein kinase 1, Cell death protein RIP, Receptor-interacting protein 1
基因 ID (Gene ID)	19766.0
蛋白 ID (SwissProt ID)	Q60855.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 75 kDa; Observed MW: 75 kDa

研究背景 (Background)

产品名称: RIP Rabbit Monoclonal Antibody
产品货号: AMRe03918




Phosphorylates DAB2IP at 'Ser-728' in a TNF-alpha-dependent manner, and thereby activates the MAP3K5-JNK apoptotic cascade (By similarity).Required for ZBP1-induced NF-kappa-B activation in response to DNA damage (PubMed:12654725, PubMed:19590578).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of RIP in RAW 264.7,NIH/3T3 lysates using RIP antibody.

注意事项 (Note)

For research use only .