

产品名称: TINF2 Rabbit Monoclonal Antibody
产品货号: AMRe03921



产品概述 (Summary)

产品名称 (Production Name)	TINF2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TINF2
别名 (Alternative Names)	DKCA3; Tin 2; TIN2; TINF 2; Tinf2
基因 ID (Gene ID)	26277
蛋白 ID (SwissProt ID)	Q9BSI4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:50 kDa;Observed MW: 50 kDa

研究背景 (Background)

Component of the shelterin complex (telosome) that is involved in the regulation of telomere length and protection.

产品名称: TINF2 Rabbit Monoclonal Antibody
产品货号: AMRe03921




Shelterin associates with arrays of double-stranded TTAGGG repeats added by telomerase and protects chromosome ends; without its protective activity, telomeres are no longer hidden from the DNA damage surveillance and chromosome ends are inappropriately processed by DNA repair pathways. Plays a role in shelterin complex assembly.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of TINF2 in MCF7, HUVEC lysates using TINF2 antibody.

注意事项 (Note)

For research use only .