

产品名称: Angiotensinogen Rabbit Monoclonal Antibody
产品货号: AMRe03924



产品概述 (Summary)

产品名称 (Production Name)	Angiotensinogen Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGT
别名 (Alternative Names)	AGT; SERPINA8; Angiotensinogen; Serpin A8
基因 ID (Gene ID)	183
蛋白 ID (SwissProt ID)	P01019.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:52 kDa;Observed MW: 52 kDa

研究背景 (Background)

Angiotensin is formed from a precursor, angiotensinogen, which is produced by the liver and found in the α -globulin

产品名称: Angiotensinogen Rabbit Monoclonal Antibody
产品货号: AMRe03924




fraction of plasma. Angiotensin stimulation through AT1 has been shown to activate the JAK/Stat pathway involving a direct interaction between JAK2 and AT1 as demonstrated by co-immunoprecipitation.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Angiotensinogen in PC-12 lysates using Angiotensinogen antibody.

注意事项 (Note)

For research use only .