

产品名称: PIM2 Rabbit Monoclonal Antibody
产品货号: AMRe03925



产品概述 (Summary)

产品名称 (Production Name)	PIM2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PIM2
别名 (Alternative Names)	PIM2; Serine/threonine-protein kinase pim-2; Pim-2h
基因 ID (Gene ID)	11040
蛋白 ID (SwissProt ID)	Q9P1W9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:34 kDa;Observed MW: 34-40 kDa

研究背景 (Background)

Promotes cell survival in response to a variety of proliferative signals via positive regulation of the I-kappaB kinase/NF-

产品名称: PIM2 Rabbit Monoclonal Antibody
产品货号: AMRe03925




kappaB cascade; this process requires phosphorylation of MAP3K8/COT. Prevents apoptosis induced by growth factor withdrawal via inhibition of caspase-3 activation, and via phosphorylation of pro-apoptotic proteins.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PIM2 in Jurkat, K562 lysates using PIM2 antibody.

注意事项 (Note)

For research use only .