

产品名称: DTN-A Rabbit Monoclonal Antibody
产品货号: AMRe03927



产品概述 (Summary)

产品名称 (Production Name)	DTN-A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DTNA
别名 (Alternative Names)	DTN; DRP3; DTN-A; Alpha-dystrobrevin; Dystrophin-related protein 3,Dystrobrevin alpha
基因 ID (Gene ID)	1837.0
蛋白 ID (SwissProt ID)	Q9Y4J8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:84 kDa;Observed MW: 84 kDa

研究背景 (Background)

产品名称: DTN-A Rabbit Monoclonal Antibody
产品货号: AMRe03927




May be involved in the formation and stability of synapses as well as being involved in the clustering of nicotinic acetylcholine receptors.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of DTN-A in Human fetal lysates using DTN-A antibody.

注意事项 (Note)

For research use only .