

产品名称: PPAR gamma Rabbit Monoclonal Antibody
产品货号: AMRe03928



产品概述 (Summary)

产品名称 (Production Name)	PPAR gamma Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPARG
别名 (Alternative Names)	CIMT1; GLM1; NR1C3; PPAR gamma; PARG; PPARG1; PPARG2;
基因 ID (Gene ID)	5468
蛋白 ID (SwissProt ID)	P37231.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 58 kDa; Observed MW: 58 kDa

研究背景 (Background)

Receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the

产品名称: PPAR gamma Rabbit Monoclonal Antibody
产品货号: AMRe03928




receptor binds to a promoter element in the gene for acyl-CoA oxidase and activates its transcription. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of PPAR gamma in Human fetal lysates using PPAR gamma antibody.

注意事项 (Note)

For research use only .