

产品名称: AP2S1 Rabbit Monoclonal Antibody
产品货号: AMRe03929



产品概述 (Summary)

产品名称 (Production Name)	AP2S1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AP2S1
别名 (Alternative Names)	AP17; AP17 delta; Ap2s1; CLAPS2; Sigma2 adaptin
基因 ID (Gene ID)	1175
蛋白 ID (SwissProt ID)	P53680.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:17 kDa;Observed MW: 17 kDa

研究背景 (Background)

Component of the adaptor protein complex 2 (AP-2). Adaptor protein complexes function in protein Transport via

产品名称: AP2S1 Rabbit Monoclonal Antibody
产品货号: AMRe03929




Transport vesicles in different membrane traffic pathways. Adaptor protein complexes are vesicle coat components and appear to be involved in cargo selection and vesicle formation.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of AP2S1 in SH-SY5Y ,293T lysates using AP2S1 antibody.

注意事项 (Note)

For research use only .