

产品名称: NSMase Rabbit Monoclonal Antibody
产品货号: AMRe03936



产品概述 (Summary)

产品名称 (Production Name)	NSMase Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SMPD2
别名 (Alternative Names)	sphingomyelin phosphodiesterase 2; ISC1; NSMASE; NSMASE1
基因 ID (Gene ID)	6610
蛋白 ID (SwissProt ID)	O60906.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa;Observed MW: 50 kDa

研究背景 (Background)

Catalyzes, at least in vitro, the hydrolysis of sphingomyelin to form ceramide and phosphocholine.


产品名称: NSMase Rabbit Monoclonal Antibody
产品货号: AMRe03936



研究领域 (Research Area)

Neuroscience; Neurology process; Neurodegenerative disease; Alzheimer's disease; Cell Type Marker; Glia marker; Oligodendrocyte marker

图片 (Image Data)



Western blot analysis of NSMase in K562 lysates using NSMase antibody.

注意事项 (Note)

For research use only .