

产品名称: Cytokeratin 14 Rabbit Monoclonal Antibody
产品货号: AMRe03940



产品概述 (Summary)

产品名称 (Production Name)	Cytokeratin 14 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KRT14
别名 (Alternative Names)	KRT14; Keratin; type I cytoskeletal 14; Cytokeratin-14; CK-14; Keratin-14; K14
基因 ID (Gene ID)	3861
蛋白 ID (SwissProt ID)	P02533.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:52 kDa;Observed MW: 52 kDa

研究背景 (Background)

The nonhelical tail domain is involved in promoting KRT5-KRT14 filaments to self-organize into large bundles and

产品名称: Cytokeratin 14 Rabbit Monoclonal Antibody
产品货号: AMRe03940




enhances the mechanical properties involved in resilience of keratin intermediate filaments in vitro.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Cytokeratin 14 in human squamous lung carcinoma, A431 lysates using Cytokeratin 14 antibody.

注意事项 (Note)

For research use only .