

产品名称: ADAR1 Rabbit Monoclonal Antibody
产品货号: AMRe03943



产品概述 (Summary)

产品名称 (Production Name)	ADAR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADAR
别名 (Alternative Names)	ADAR; Adar1; AGS6; DRADA; Dsh; Dsrad; IFI4; P136
基因 ID (Gene ID)	103
蛋白 ID (SwissProt ID)	P55265.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:136 kDa;Observed MW: 150 kDa

研究背景 (Background)

Converts multiple adenosines to inosines and creates I/U mismatched base pairs in double-helical RNA substrates without

产品名称: ADAR1 Rabbit Monoclonal Antibody
产品货号: AMRe03943




apparent sequence specificity.

研究领域 (Research Area)

Microbiology

图片 (Image Data)



Western blot analysis of ADAR1 in HeLa lysates using ADAR1 antibody.

注意事项 (Note)

For research use only .