

产品名称: PGC1 beta Rabbit Monoclonal Antibody
产品货号: AMRe03945



产品概述 (Summary)

产品名称 (Production Name)	PGC1 beta Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPARGC1B
别名 (Alternative Names)	PERC; PGC1; PPARGC1; Ppargc1b
基因 ID (Gene ID)	133522
蛋白 ID (SwissProt ID)	Q86YN6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:113 kDa;Observed MW: 113 kDa

研究背景 (Background)

Plays a role of stimulator of transcription factors and nuclear receptors activities. Activates transcriptional activity of estrogen

产品名称: PGC1 beta Rabbit Monoclonal Antibody
产品货号: AMRe03945




receptor alpha, nuclear respiratory factor 1 (NRF1) and glucocorticoid receptor in the presence of glucocorticoids.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PGC1 beta in HeLa, MCF-7 lysates using PGC1 beta antibody.

注意事项 (Note)

For research use only .