

产品名称: Emi1 Rabbit Monoclonal Antibody
产品货号: AMRe03948



产品概述 (Summary)

产品名称 (Production Name)	Emi1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FBXO5
别名 (Alternative Names)	EMI1; FBX5; Fbxo31; fbxo5
基因 ID (Gene ID)	26271
蛋白 ID (SwissProt ID)	Q9UKT4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:50 kDa;Observed MW: 56kDa

研究背景 (Background)

Regulates progression through early mitosis by inhibiting the anaphase promoting complex/cyclosome (APC). Binds to the

产品名称: Emi1 Rabbit Monoclonal Antibody
产品货号: AMRe03948




APC activators CDC20 and FZR1/CDH1 to prevent APC activation. Can also bind directly to the APC to inhibit substrate-binding.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Emi1 in HeLa, 293 lysates using Emi1 antibody.

注意事项 (Note)

For research use only .