

产品名称: PPIC Rabbit Monoclonal Antibody
产品货号: AMRe03949



产品概述 (Summary)

产品名称 (Production Name)	PPIC Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPIC
别名 (Alternative Names)	PPIC,CYPC,Rotamase C,Cyclophilin C,Peptidyl-prolyl cis-trans isomerase C
基因 ID (Gene ID)	5480.0
蛋白 ID (SwissProt ID)	P45877.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa;Observed MW: 23 kDa

研究背景 (Background)

PPIase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist

产品名称: PPIC Rabbit Monoclonal Antibody
产品货号: AMRe03949




protein folding.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of PPIC in HCT116 lysates using PPIC antibody.

注意事项 (Note)

For research use only .