

产品名称: CAD Rabbit Monoclonal Antibody
产品货号: AMRe03951



产品概述 (Summary)

产品名称 (Production Name)	CAD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAD
别名 (Alternative Names)	CAD; CAD protein
基因 ID (Gene ID)	790
蛋白 ID (SwissProt ID)	P27708.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:243 kDa;Observed MW: 243 kDa

研究背景 (Background)

Carbamoyl phosphate synthetase-aspartate carbamoyltransferase-dihydroorotase (CAD) is a multifunctional protein that

产品名称: CAD Rabbit Monoclonal Antibody
产品货号: AMRe03951




initiates and regulates mammalian de novo pyrimidine biosynthesis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of CAD in 293 lysates using CAD antibody.

注意事项 (Note)

For research use only .