

产品名称: Gephyrin Rabbit Monoclonal Antibody
产品货号: AMRe03956



产品概述 (Summary)

产品名称 (Production Name)	Gephyrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPHN
别名 (Alternative Names)	Domain E; Domain G; GEPH; GPH; GPHN; GPHRYN
基因 ID (Gene ID)	10243
蛋白 ID (SwissProt ID)	Q9NQX3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:80 kDa;Observed MW: 93 kDa

研究背景 (Background)

Microtubule-associated protein involved in membrane protein-cytoskeleton interactions. It is thought to anchor the

产品名称: Gephyrin Rabbit Monoclonal Antibody
产品货号: AMRe03956




inhibitory glycine receptor (GLYR) to subsynaptic microtubules (By similarity). Catalyzes two steps in the biosynthesis of the molybdenum cofactor.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Gephyrin in 293T lysates using Gephyrin antibody.

注意事项 (Note)

For research use only .