

产品名称: GRB7 Rabbit Monoclonal Antibody
产品货号: AMRe03957



产品概述 (Summary)

产品名称 (Production Name)	GRB7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRB7
别名 (Alternative Names)	B47; Epidermal growth factor receptor GRB 7; GRB7 adapter protein; GRB7
基因 ID (Gene ID)	2886
蛋白 ID (SwissProt ID)	Q14451.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:60 kDa;Observed MW: 60 kDa

研究背景 (Background)

Interacts with the cytoplasmic domain of the epidermal growth factor receptor which is then inhibited. The interaction is

产品名称: GRB7 Rabbit Monoclonal Antibody
产品货号: AMRe03957




mediated by the SH2 domain. Also binds to ERBB2.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GRB7 in HepG2 lysates using GRB7 antibody.

注意事项 (Note)

For research use only .