

产品名称: NIPP1 Rabbit Monoclonal Antibody
产品货号: AMRe03966



产品概述 (Summary)

产品名称 (Production Name)	NIPP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP1R8
别名 (Alternative Names)	ARD1; ARD-1; NIPP1; NIPP-1; PRO2047
基因 ID (Gene ID)	5511.0
蛋白 ID (SwissProt ID)	Q12972.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:39 kDa;Observed MW: 39 kDa

研究背景 (Background)

Inhibitor subunit of the major nuclear protein phosphatase-1 (PP-1). It has RNA-binding activity but does not cleave RNA

产品名称: NIPP1 Rabbit Monoclonal Antibody
产品货号: AMRe03966




and may target PP-1 to RNA-associated substrates. May also be involved in pre-mRNA splicing. Binds DNA and might act as a transcriptional repressor. Seems to be required for cell proliferation.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of NIPP1 in HeLa, Jurkat lysates using NIPP1 antibody.

注意事项 (Note)

For research use only .