

产品名称: GARNL1 Rabbit Monoclonal Antibody
产品货号: AMRe03975



产品概述 (Summary)

产品名称 (Production Name)	GARNL1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RALGAPA1
别名 (Alternative Names)	p240; GRIPE; GARNL1; TULIP1; RalGAPalpha1
基因 ID (Gene ID)	253959
蛋白 ID (SwissProt ID)	Q6GYQ0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:230 kDa;Observed MW: 230 kDa

研究背景 (Background)

Catalytic subunit of the heterodimeric RalGAP1 complex which acts as a GTPase activator for the Ras-like small GTPases

产品名称: GARNL1 Rabbit Monoclonal Antibody
产品货号: AMRe03975




RALA and RALB.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of GARNL1 in 293, lysates using GARNL1 antibody.

注意事项 (Note)

For research use only .