

产品名称: RFC2 Rabbit Monoclonal Antibody
产品货号: AMRe03977



产品概述 (Summary)

产品名称 (Production Name)	RFC2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% protective protein.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RFC2 RFC2, Replication factor C subunit 2, Activator 1 40 kDa subunit (A1 40 kDa subunit), Activator 1 subunit 2, Replication factor C 40 kDa subunit (RF-C 40 kDa subunit; RFC40)
别名 (Alternative Names)	
基因 ID (Gene ID)	5982.0
蛋白 ID (SwissProt ID)	P35250.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

产品名称: RFC2 Rabbit Monoclonal Antibody
产品货号: AMRe03977




研究背景 (Background)

The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins proliferating cell nuclear antigen (PCNA) and activator 1. This subunit binds ATP (By similarity).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of RFC2 in K562, HeLa lysates using RFC2 antibody.

注意事项 (Note)

For research use only .